

NEO DELTA SELF SAFE 1 NEO DELTA SELF SAFE T

STERILE EO

EN - INTENDED USE.....	3
IT - DESTINAZIONE D'USO.....	5
FR - UTILISATION PRÉVUE.....	7
DE - GEBRAUCHSBESTIMMUNG	9
ES - USO PREVISTO.....	11
CS - ZAMÝŠLENÉ POUŽITÍ.....	13
DA - TILSIGTET BRUG	15
FI - KÄYTTÖTARKOITUS.....	17
EL - ΠΡΟΒΛΕΠΟΜΕΝΗ ΧΡΗΣΗ.....	19
HU - RENDELTETÉSSZERŰ HASZNÁLAT	21
NL - BEOOGD GEBRUIK	23
NO - TILTENKT BRUK	25
PL - PRZEZNACZENIE	27
PT - USO PREVISTO.....	29
RO - UTILIZARE PREVĂZUTĂ	31
SK - URČENÉ POUŽITIE	33
LV - PAREDZĒTAIS LIETOJUMS.....	35
LT - NUMATYTOJI PASKIRTIS	37
BG - ПРЕДНАЗНАЧЕНИЕ.....	39
SI - PREDVIDENA UPORABA.....	41
TR - KULLANIM AMACI.....	43
RU - НАЗНАЧЕНИЕ	45

EN - INTENDED USE

Peripheral Venous catheter provided with passive safety mechanism. Peripheral IV Catheter for short term peripheral venous access that, in combination with other medical device, allows the administration of fluids. Once placed in the vein, the IV Catheter can be connected to others medical devices to administer therapeutical solutions or drugs. IV Catheters can be connected with others medical devices provided with luer lock or luer fitting connections such as infusion sets, extension lines, stopcocks and syringes. In case of high pressure administration of fluids, like contrast media, the device is connected to extension lines connected to power injectors to administer the contrast media at high pressure.

The Catheter can be used on any patient population while taking into account the vascular anatomy of the patient and of the adequacy of the procedure.

The catheters are suitable for use with pressure injectors (max. 325psi)

DO NOT USE 26G IV CATHETERS FOR HIGH PRESSURE TREATMENTS.

CONTRAINDICATION

Do not use on patients with known hypersensitivity to the materials used.

MATERIALS USED

Stainless Steel, radiopaque Polyurethane (PUR), Polypropylene (PP), Polyoxymethylene (POM). Lubricants: polymethylsiloxane and aminofunctional polymethylsiloxane

GENERAL WARNINGS

- Choose the correct IV Cannula size.
- Open the primary packaging following the indications (PEEL).
- Intended for single patient use only.
- Do not resterilize.
- DO NOT RE-USE. Re-use may lead to infection or other illness/injury and may compromise device functionality.
- Exposure to blood, either through percutaneous puncture with a contaminated needle or via mucous membranes or nonintact skin, may lead to serious illness such as hepatitis, HIV (AIDS), or other infectious diseases.
- Do not attempt to override or defeat the needle shielding mechanism as a contaminated needlestick injury may occur.
- Do not leave Luer slip connections unattended as they may become disconnected and result in blood exposure or blood loss.
- After use, discard the device according to your facility policy.
- Use the device to administer proper solutions and fluids in order to reduce as far as possible phlebitis caused by hypotonic or hypertonic solution.
- IV catheters are intended to be used only by trained and qualified personnel for using this kind of medical device and for their management

GENERAL PRECAUTIONS

- Use these catheters in conjunction with your facility policy.
- Ensure aseptic technique, proper skin preparation, and continued protection of the site consistent with accepted standards of practice and your facility policy.
- Replace catheter according to your facility policy, relevant guidelines, or if the integrity of the device has been compromised.
- Follow Universal/Standard Precautions when inserting, maintaining, and removing these catheters to prevent exposure to bloodborne pathogens.
- Adhere to all the contraindications, warnings, precautions, and instruction for all infusate, as specified by their manufacturer
- Do not bend the needle before or during use.
- Do not use scissors or other sharp instruments at or near the insertion site.
- Use only with ISO 80369-7 compliant Luer connections. Non-ISO compliant Luer connections may cause leakage.
- Do not overtighten ISO 80369-7 Luer connections as damage may occur.

INSTRUCTIONS


1. **CAUTION:** Do not use if package is damaged or opened or if the expiration date has passed.
2. Prepare insertion site according to your facility policy.
3. Open package and remove catheter by grasping the grip.
4. Holding the wings, twist and remove needle cover in a straight, outward motion. Discard needle cover according to your facility policy.
5. Inspect device for damage before insertion. **CAUTION:** Do not use if device is damaged.
6. Access vessel using a low angle of insertion. **WARNING:** If the needle is partially or completely withdrawn from the catheter tubing during insertion, do not re-insert the needle into the catheter tubing as damage may occur.
7. Observe blood return in the flash chamber
8. Advance the catheter off the needle into the vessel. **CAUTION:** Do not pull back on the needle during catheter advancement.
9. Stabilize the catheter by gently pressing the wings to the skin and apply venous compression. **WARNING:** Venous compression is required to prevent blood leakage from the catheter hub.

10. Pull the needle out of the catheter until an audible “click” is heard indicating that the needle safety shield has activated.
WARNING: If the needle fails to shield, immediately dispose of the needle in a puncture resistant, leak-proof sharps container keeping the needle point away from the body and finger at all times.
11. Continue to apply venous compression and detach shielded needle from the catheter hub.
12. Remove Luer lock cap before discarding shielded needle into a puncture resistant, leak-proof sharps container.
13. Securely connect the Luer lock cap or other Luer device to the catheter hub.
14. Stabilize catheter and apply a sterile dressing according to your facility policy.
15. Flush or begin infusion.
16. Perform vascular access site care and dressing changes according to your facility policy.
17. Upon removal, examine catheter to ensure it is intact and discard according to your facility policy

Warnings. (PRECAUTIONS.)

- Read the instructions before use.
- Use protective gloves.
- The product must be used immediately after the packaging has been opened.
- In case of incorrect transport and/or manipulation, the device or packaging could be subject to structural and/or functional damage.
- Do not use scissors at or near the insertion site.
- **If the device is used at a high pressure or with injectors:**
 - Connect directly to the pressure infusion system with the end luer-lock connector of the device.
 - Remove all the accessories connected to the device and replace them with a luer-lock cap where necessary.
 - Always check the patency of the device before use.
 - Never exceed the maximum pressure of 325 psi
- Flush the device immediately after the administration of medicines or biological fluids.
- Immediately remove any needle that has no coating, always keeping the tip away from your body and fingers.
- Do not expose to heat or direct sunlight.

Note: Any serious incident that has occurred in relation to the device should be reported to the Manufacturer and the Competent Authority of the Member State in which the user and /or the patient is established.

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