

Saflo Microvolume Extension Set

Infutech Ltd are pleased to present the first of the new Saflo™range of products designed specifically for use in a wide range of subcutaneous infusion applications.

The products in this range have been created following extensive work with practitioners in this area to meet their daily requirements. The Saflo products now extend to a selection of different needle safe infusion set configurations available directly from us or through our distributors.

Code	Box Qty	Description
SFX100-AS	10	Saflo micro-volume extension set 100 cm with anti-siphon valve
SFX150-00	10	Saflo micro-volume extension set 150 cm
SFX-102Y	10	Micro-volume extension set Y configuration 100 cm



Standard set 150cm



Anti siphon set 100cm



Y-Set 100cm



The principle features of this new range of extension sets are:

- Very low priming volumes: approximately 0.20mL for standard and Y-Set approximately 0.15mL for anti siphon set
- Tri-layer tubing construction using DEHP free PVC for a soft outer layer for comfort with a polyethylene inner layer for maximum drug compatibility
- Kink resistant tubing
- Compatible with all standard infusion pumps