

Saflo Product Disposal and Recycling Information.

Please consult your local authority's guidance on recycling in your local area as this may differ from the information provided below.

Once the product has been used it must be treated as contaminated waste and disposed of according to the local hospital instructions.

The following applies to the packaging materials for the products.

Saflo 90 Catheter Products Codes: SF90-109, SF90-106, SF90-PO9, SF90-P06

Saflo Catheter Product Code SFC-001

Item	Material	Recycling Information
Product Blister Pack	PET G	PETG is fully recyclable and can be disposed of in regular plastic recycling bins
Blister Pack Lid	TYVEK/Paper	Tyvek is manufactured from a High Density Polyethelene (HDPE) material and is fully recyclable. Dupont, the manufacturers of the TYVEK brand, have established a UK recycling centre at: Spenic Converters LLP Unit 11 Caker Stream Road Mill Lane, Alton Hampshire GU34 2QA www.spenic-converters.co.uk
Product IFU	Paper	Disposal in regular disposal facilities
Delivery Packaging	Card/Paper	Disposal in regular disposal facilities
Shipping Cartons	Cardboard	Disposal in regular disposal facilities

Saflo Tubing Codes: SFT100-00, SFT060.00, SFT100-SP, SFP-001

Item	Material	Recycling Information
Product Pack Clear Material	PE	PE is fully recyclable and can be disposed of in regular plastic recycling bins
Product Pack Paper	Paper	Disposal in regular disposal facilities
Product IFU	Paper	Disposal in regular disposal facilities
Delivery Packaging	Card/Paper	Disposal in regular disposal facilities
Shipping Cartons	Cardboard	Disposal in regular disposal facilities

Please contact us if you have any particular concerns concerning disposal.