

DELTA SELF SAFE

Safety Peripheral Intravenous Cather

DELTA SELFSAFE® **FEP**

DELTA SELFSAFE® **PUR**



 **DELTA MED**
MEDICAL DEVICES
Solutions for Health

PART OF



NUMANTEC
HUMAN INNOVATION TECHNOLOGY

PERIPHERAL INTRAVENOUS CATHETERS

SAFETY PERIPHERAL INTRAVENOUS CATHETER

DELTA SELFSAFE® **FEP** DELTA SELFSAFE® **PUR**

- Catheter material: **clear**; flexible; **anti kinking**; **visible** under X-Ray (3 radiopaque lines).
- Available in **Polyurethane** (PUR) version (particularly heat-sensitive: softening once in vein) and **Teflon** (FEP) version.
- Back-cut, sharp needle bevel for smooth, **non-traumatic insertion**.
- Color coded catheter hub/port for **immediate identification** of catheter **gauges**.
- Ultra clear **Flashback Chamber** for **quick visualization of blood return** and confirmation of **vein access**.



PASSIVE SAFETY DEVICE

- It automatically activates while withdrawing the needle
- Audible **CLICK** confirming successful activation
- Double protection mechanism: inside clip covering the needle tip plus outside housing protecting further from scratches and blood exposure

Neo Delta Self Safe 2 / Delta Self Safe 2

PORTED AND WINGED CATHETER



Neo DELTA SELF SAFE 2 (PUR)	DELTA SELF SAFE 2 (FEP)	DESCRIPTION	GAUGE	CATH LENGTH	CATH. OD	CATH. ID	FLOW RATE (ml/min)	NEEDLE OD	SALES UNIT / UNIT PER BOX
3704022	-	Ported and winged 26G	26	19 mm (3/4")	0.64 mm	0.45 mm	19	0.4 mm	50/200
3713122	3711522	Ported and winged 24G	24	19 mm (3/4")	0.74 mm	0.55 mm	29	0.5 mm	50/200
3723122	3721522	Ported and winged 22G	22	25 mm (1")	0.90 mm	0.65 mm	42	0.6 mm	50/200
3733122	3731522	Ported and winged 20G	20	32 mm (1" 1/4)	1.00 mm	0.75 mm	59	0.7 mm	50/200
3783122	3781522	Ported and winged 18G	18	32 mm (1" 1/4)	1.30 mm	0.95 mm	103	0.9 mm	50/200
3743122	3741522	Ported and winged 18G	18	45 mm (1" 3/4)	1.30 mm	0.95 mm	96	0.9 mm	50/200
3753122	-	Ported and winged 17G	17	45 mm (1" 3/4)	1.50 mm	1.15 mm	155	1.1 mm	50/200
3763122	-	Ported and winged 16G	16	45 mm (1" 3/4)	1.75 mm	1.35 mm	225	1.3 mm	50/200
3773122	-	Ported and winged 14G	14	45 mm (1" 3/4)	2.00 mm	1.55 mm	290	1.5 mm	50/200

Neo Delta Self Safe 1

STRAIGHT WINGED CATHETER



Neo DELTA SELF SAFE 1 (PUR)	DESCRIPTION	GAUGE	CATH LENGTH	CATH. OD	CATH. ID	FLOW RATE (ml/min)	NEEDLE OD	SALES UNIT / UNIT PER BOX
3708122	Straight winged 26G	26	19 mm (3/4")	0.64 mm	0.45 mm	19	0.4 mm	50/200
3718122	Straight winged 24G	24	19 mm (3/4")	0.74 mm	0.55 mm	29	0.5 mm	50/200
3728122	Straight winged 22G	22	25 mm (1")	0.90 mm	0.65 mm	42	0.6 mm	50/200
3738122	Straight winged 20G	20	32 mm (1" 1/4)	1.00 mm	0.75 mm	59	0.7 mm	50/200
3788122	Straight winged 18G	18	32 mm (1" 1/4)	1.30 mm	0.95 mm	103	0.9 mm	50/200
3748122	Straight winged 18G	18	45 mm (1" 3/4)	1.30 mm	0.95 mm	96	0.9 mm	50/200
3768122	Straight winged 16G	16	45 mm (1" 3/4)	1.75 mm	1.35 mm	225	1.3 mm	50/200
3778122	Straight winged 14G	14	45 mm (1" 3/4)	2.00 mm	1.55 mm	290	1.5 mm	50/200

Neo Delta Self Safe T

STRAIGHT CATHETER



Neo DELTA SELF SAFET (PUR)	DESCRIPTION	GAUGE	CATH LENGTH	CATH. OD	CATH. ID	FLOW RATE (ml/min)	NEEDLE OD	SALES UNIT / UNIT PER BOX
3708522	Straight 26G	26	19 mm (3/4")	0.64 mm	0.45 mm	19	0.4 mm	50/200
3718522	Straight 24G	24	19 mm (3/4")	0.74 mm	0.55 mm	29	0.5 mm	50/200
3728522	Straight 22G	22	25 mm (1")	0.90 mm	0.65 mm	42	0.6 mm	50/200
3738522	Straight 20G	20	32 mm (1" 1/4)	1.00 mm	0.75 mm	59	0.7 mm	50/200
3788522	Straight 18G	18	32 mm (1" 1/4)	1.30 mm	0.95 mm	103	0.9 mm	50/200
3748522	Straight 18G	18	45 mm (1" 3/4)	1.30 mm	0.95 mm	96	0.9 mm	50/200
3768522	Straight 16G	16	45 mm (1" 3/4)	1.75 mm	1.35 mm	225	1.3 mm	50/200
3778522	Straight 14G	14	45 mm (1" 3/4)	2.00 mm	1.55 mm	290	1.5 mm	50/200