

DELTA VEN

Closed System Safety Peripheral Intravenous Catheter





PERIPHERAL INTRAVENOUS CATHETERS

CLOSED SYSTEM SAFETY PIV C FOR VENOUS AND SUBCUTANEOUS ACCESS

The Deltaven is the **only 26G neonatal closed system IV catheter**, designed to meet the unique needs of the smallest veins.



Flexible stabilization platform with gentle grooves enhancing securement and increasing patient's comfort



*The second needle free connector will be capped and loose in the unitary blister



CODE	DESCRIPTION	GAUGE	CATHETER O.D. (mm)	CATHETER LENGTH (mm)	GRAVITY FLOW RATE (ml/min)	MAXIMUM POWER INJECTOR FLOW RATE (ml/sec)*	SALES UNIT / UNIT PER BOX
3904722	Single Port	26	0.64	19	15	Not for use with power injectors	25/100
3904822	Single Port w/stopcock and needle-free connector	26	0.64	19	15		20/80
3904922	Single Port w/stopcock	26	0.64	19	15		20/80
3914722	Single Port	24	0.74	19	24	4	25/100
3914822	Single Port w/stopcock and needle-free connector	24	0.74	19	24	4	20/80
3914922	Single Port w/stopcock	24	0.74	19	24	4	20/80
3924722	Single Port	22	0.90	19	37	5	25/100
3924822	Single Port w/stopcock and needle-free connector	22	0.90	19	37	5	20/80
3924922	Single Port w/stopcock	22	0.90	19	37	5	20/80
3925722	Single Port	22	0.90	25	31	5	25/100
3925822	Single Port w/stopcock and needle-free connector	22	0.90	25	31	5	20/80
3925922	Single Port w/stopcock	22	0.90	25	31	5	20/80
3935722	Single Port	20	1.00	25	59	8	25/100
3935822	Single Port w/stopcock and needle-free connector	20	1.00	25	59	8	20/80
3935922	Single Port w/stopcock	20	1.00	25	59	8	20/80
3936722	Single Port	20	1.00	32	56	8	25/100
3936822	Single Port w/stopcock and needle-free connector	20	1.00	32	56	8	20/80
3936922	Single Port w/stopcock	20	1.00	32	56	8	20/80
3937722	Single Port	20	1.00	45	51	7	25/100
3937822	Single Port w/stopcock and needle-free connector	20	1.00	45	51	7	20/80
3937922	Single Port w/stopcock	20	1.00	45	51	7	20/80
3945722	Single Port	18	1.30	25	86	15	25/100
3945822	Single Port w/stopcock and needle-free connector	18	1.30	25	86	15	20/80
3945922	Single Port w/stopcock	18	1.30	25	86	15	20/80
3946722	Single Port	18	1.30	32	83	14	25/100
3946822	Single Port w/stopcock and needle-free connector	18	1.30	32	83	14	20/80
3946922	Single Port w/stopcock	18	1.30	32	83	14	20/80
3947722	Single Port	18	1.30	45	78	13	25/100
3947822	Single Port w/stopcock and needle-free connector	18	1.30	45	78	13	20/80
3947922	Single Port w/stopcock	18	1.30	45	78	13	20/80
3956722	Single Port	16	1.75	32	125	20	25/100
3956822	Single Port w/stopcock and needle-free connector	16	1.75	32	125	20	20/80
3956922	Single Port w/stopcock	16	1.75	32	125	20	20/80

*Maximum Power Injector Pressure Limit Setting: 325 psi. Contrast media used: IOMERON 400 featuring 27.5 mPa.s viscosity at 20°C

CODE	DESCRIPTION	GAUGE	CATHETER O.D. (mm)	CATHETER LENGTH (mm)	GRAVITY FLOW RATE (ml/min)	MAXIMUM POWER INJECTOR FLOW RATE (ml/sec)*	SALES UNIT / UNIT PER BOX
3900522	Dual Port	26	0.64	19	15	Not for use with power injectors	25/100
3900622	Dual Port w/needle-free connector	26	0.64	19	15		25/100
3900422	Dual Port w/dual needle-free connector	26	0.64	19	15		20/80
3911022	Dual Port	24	0.74	19	24	4	25/100
3910622	Dual Port w/needle-free connector	24	0.74	19	24	4	25/100
3911422	Dual Port w/dual needle-free connector	24	0.74	19	24	4	20/80
3920522	Dual Port	22	0.90	19	37	5	25/100
3920622	Dual Port w/needle-free connector	22	0.90	19	37	5	25/100
3920422	Dual Port w/dual needle-free connector	22	0.90	19	37	5	20/80
3921522	Dual Port	22	0.90	25	31	5	25/100
3921622	Dual Port w/needle-free connector	22	0.90	25	31	5	25/100
3921422	Dual Port w/dual needle-free connector	22	0.90	25	31	5	20/80
3931522	Dual Port	20	1.00	25	59	8	25/100
3931622	Dual Port w/needle-free connector	20	1.00	25	59	8	25/100
3931422	Dual Port w/dual needle-free connector	20	1.00	25	59	8	20/80
3932522	Dual Port	20	1.00	32	56	8	25/100
3932622	Dual Port w/needle-free connector	20	1.00	32	56	8	25/100
3932422	Dual Port w/dual needle-free connector	20	1.00	32	56	8	20/80
3933522	Dual Port	20	1.00	45	51	7	25/100
3933622	Dual Port w/needle-free connector	20	1.00	45	51	7	25/100
3933422	Dual Port w/dual needle-free connector	20	1.00	45	51	7	20/80
3941522	Dual Port	18	1.30	25	86	15	25/100
3941622	Dual Port w/needle-free connector	18	1.30	25	86	15	25/100
3941422	Dual Port w/dual needle-free connector	18	1.30	25	86	15	20/80
3942522	Dual Port	18	1.30	32	83	14	25/100
3942622	Dual Port w/needle-free connector	18	1.30	32	83	14	25/100
3942422	Dual Port w/dual needle-free connector	18	1.30	32	83	14	20/80
3943522	Dual Port	18	1.30	45	78	13	25/100
3943622	Dual Port w/needle-free connector	18	1.30	45	78	13	25/100
3943422	Dual Port w/dual needle-free connector	18	1.30	45	78	13	20/80
3952522	Dual Port	16	1.75	32	125	20	25/100
3952622	Dual Port w/needle-free connector	16	1.75	32	125	20	25/100
3952422	Dual Port w/dual needle-free connector	16	1.75	32	125	20	20/80

*Maximum Power Injector Pressure Limit Setting: 325 psi. Contrast media used: IOMERON 400 featuring 27.5 mPa.s viscosity at 20°C